

▼ **PULL** ▼

HEARTSTART
DEFIBRILLATORS



55+ lbs / 25+ kg

HeartStart First Aid Defibrillator

(This image represents the actual size of the HeartStart First Aid Defibrillator)

**HEART
START**
HOME DEFIBRILLATOR

PHILIPS

HEARTSTART

DEFIBRILLATORS

HeartStart First Aid Defibrillator

Each year sudden cardiac arrest strikes approximately 30,000¹ people in Australia & New Zealand. For every minute defibrillation is delayed the chance of survival decreases by about 10%. Defibrillation is part of the Australian and New Zealand Resuscitation Council Guidelines.

The HeartStart First Aid is designed for use by everyone, at any time, at any location. It is fast to learn, easy to use and always ready. **Ideal for the office, sporting clubs, shopping centres, schools and many other public places.**



- **Fast to Learn:** specifically designed for the lay responder.^{2,3}
- **Easy to Use.**^{2,3}
- **Always Ready** with daily self testing including: system electronics, battery and pads readiness.
- **CPR assistance** for both Adults and Children with clear Australasian Voice Prompts.
- **Quick Shock Capability.**
- **Long Life Battery** that is easy to replace (200 shocks, 4 year standby and 5 year shelf life).
- **Preconnected SMART pads.**

Quick Shock the HeartStart Advantage.

Quick Shock is a unique feature offered only by Laerdal's HeartStart Defibrillator Range. The HeartStart Defibrillator is able to deliver a shock 10 seconds or less after the end of a CPR pause.

No other automated external defibrillator is able to do this as quickly.

With more than **600,000** units deployed worldwide, more workplaces, organisations and communities place their trust in Philips/Laerdal HeartStart AEDs than any other brand when establishing early defibrillation programs. In Australia and New Zealand this includes Ambulance services, Airlines (QANTAS & Air New Zealand), Life Saving organisations (Royal Life Saving & Surf Life Saving Australia), major airports, schools and corporations.

Easy To Use: Simple as 1, 2, 3



1-Turn ON, follow voice prompts



2-Apply Pads



3-Press to Shock

References:

1. Estimate as provided by Australian Resuscitation Council.
2. According to an independent study, "across all measures of performance, behaviour and subjective experience the Philips/Laerdal AED (HeartStart First Aid) showed superior performance over the three other devices (other devices were, Medtronic (Physio Control) CR+, Zoll AED + & Cardiac Science Powerheart"). Andre AD, When Every Minute Counts: All Automatic External Defibrillators are not created equal, 2003 p.1-5.
3. Only defibrillator that is currently available OTC (over the counter) in the US. Other models only available via a prescription from a General Practitioner. Source US FDA.

FitnessFirst



MFB



Royal Life Saving
THE ROYAL LIFE SAVING SOCIETY OF AUSTRALIA

HEARTSTART

DEFIBRILLATORS

Selection of Accessories

With our range of defibrillators we also have a selection of accessories available including defibrillator training carry cases, wall brackets and user DVD's. For more technical information on our defibrillators and details on our accessories please visit www.heartstart.com.au or www.heartstart.co.nz



HeartStart First Aid Trainer with Little Anne Manikin - Trainers are available for all HeartStart Models (except MRx). M5086A



Infant/Child Pads to suit HeartStart First Aid (Pads and Infant/Child Keys also available for FRx and FR2+). M5072A



Standard Carry Case to suit HeartStart First Aid (All HeartStart First Aids sold come free with a Slimline version). 941340



Pelican Case (Waterproof) to suit HeartStart First Aid and FRx. YC




Adult Training Pads to suit HeartStart First Aid. M5073A



HeartStart Wall Box with Alarm. Battery operated alarm with sound and light when door is opened. 941343



Basic AED Cabinet. Battery operated alarm with sound when door is opened. 9898 031 36531



HeartStart Wall Bracket to suit HeartStart First Aid, FRx and FR2+. M3857A

Australia:

LAERDAL Pty. Ltd. PO Box 52 Oakleigh VIC 3166 Free Call: 1800-DEFIBS (1800 333 427) Free Fax: 1800 635 835 Email: customerservice@laerdal.com.au www.heartstart.com.au

New Zealand:

For further information on HeartStart Defibrillators or proposals for your organisation, please contact St John on: 0800 FIRST AID (0800 347 782) Email: enquiries@stjohn.org.nz or visit www.stjohn.org.nz/aed

™ Trademark of Laerdal Pty. Ltd.

Your Local Distributor:

