

HEARTSTART
DEFIBRILLATORS

ON/OFF

1



HeartStart FRx



2



PHILIPS

3



(This image represents the actual size of the HeartStart FRx Defibrillator)

HEARTSTART

DEFIBRILLATORS

HeartStart FRx Defibrillator

Each year sudden cardiac arrest strikes approximately 30,000¹ people in Australia & New Zealand. For every minute defibrillation is delayed the chance of survival decreases by about 10%. Defibrillation is part of the Australian and New Zealand Resuscitation Council Guidelines.

The HeartStart FRx is an easy to use, rugged and reliable AED for those who get there first. **With its International Protection (IP) rating of IP 55 it's the perfect AED for industrial sites, use near water and extreme environments.**



- **Specifically designed for its purpose.**
- **Easy to Use.**
- **Always Ready with Daily Self Testing including: system electronics, battery, and pads readiness.**
- **CPR assistance for both Adults and Children with clear Australasian Voice Prompts.**
- **Quick Shock Capability.**
- **Long Life Battery that is easy to replace (200 shocks, 4 year standby and 5 year shelf life).**
- **IP55 Rating which protects the unit from dust and jetting water.**
- **Preconnected SMART pads for adults and children*.**

* paediatric key not included please order separately

Quick Shock the HeartStart Advantage.

Quick Shock is a unique feature offered only by Laerdal's HeartStart Defibrillator Range. The HeartStart Defibrillator is able to deliver a shock 10 seconds or less after the end of a CPR pause.

No other automated external defibrillator is able to do this as quickly.

With more than 600,000 units deployed worldwide, more workplaces, organisations and communities place their trust in Philips/Laerdal HeartStart AEDs than any other brand when establishing early defibrillation programs. In Australia and New Zealand this includes Ambulance services, Airlines (QANTAS & Air New Zealand), Life Saving organisations (Royal Life Saving & Surf Life Saving Australia), major airports, schools and corporations.

Easy To Use: Simple as 1, 2, 3



1-Turn ON, follow voice prompts



2-Apply Pads



3-Press to Shock

References:

1. Estimate as provided by Australian Resuscitation Council.

FitnessFirst



MFB



Royal Life Saving
THE ROYAL LIFE SAVING SOCIETY AUSTRALIA

HEARTSTART

DEFIBRILLATORS

Selection of Accessories

With our range of defibrillators we also have a selection of accessories available including defibrillator trainers, carry cases, wall brackets and user DVD's. For more technical information on our defibrillators and details on our accessories please visit www.heartstart.com.au or www.heartstart.co.nz



HeartStart FRx Trainer (Trainers are available for all HeartStart Models, except MRx).
861307



Infant/Child Key to suit HeartStart FRx (Infant/Child pads also available for HeartStart First Aid and FR2+).
9898 031 39311



Carry Case to suit HeartStart FRx (Every HeartStart FRx sold comes with this case as part of the package).
941360



Pelican Case (Waterproof) to suit HeartStart First Aid and FRx.
YC



Adult Training pads to suit HeartStart FRx.
9898 031 39271



HeartStart Wall Box with Alarm. Battery operated alarm with sound and light when door is opened.
941343



Basic AED Cabinet. Battery operated alarm with sound when door is opened.
9898 031 36531



HeartStart Wall Bracket to suit HeartStart First Aid, FRx and FR2+.
M3857A

Australia:

LAERDAL Pty. Ltd. PO Box 52 Oakleigh VIC 3166 Free Call: 1800-DEFIBS (1800 333 427) Free Fax: 1800 635 835 Email: customer.service@laerdal.com.au www.heartstart.com.au

New Zealand:

For further information on HeartStart Defibrillators or proposals for your organisation, please contact St John on: 0800 FIRST AID (0800 347 782) Email: enquiries@stjohn.org.nz or visit www.stjohn.org.nz/aed

™ Trademark of Laerdal Pty. Ltd.

Your Local
Distributor: