

Defib Case

Specifications



Product Specifications

- Waterproof and dust-proof
- Watertight seal around lid
- Automatic pressure release valve
- Easy to open double throw latches
- Load tested soft grip handle
- Holes for double padlocking
- Holes for shoulder strap
- Stackable interlocking system
- Extra fixing points
- Screen printing space and rectangular resin label space on the lid
- Cubed foams in the bottom and convoluted foams in the lid available as standard
- Thick body material
- Body and handle: polypropylene copolymer
- Fully length hinges
- Nylon pins, corrosion free
- Latches: nylon
- Pins: nylon + 30% fiberglass
- Colours: Orange and black
- Weight without foam: 2760g
- Weight with foam: 3040g
- Moulding: injection molded
- Range temperature: -30/90°C (-22/194°F)
- Features: IP67

Dimensions

External dimensions: 464 x 366 x H176mm (18.27 x 14.41 x H6.93in)

Internal dimensions: 426 x 290 x H159mm (16.77 x 11.42 x H6.26in)

Lid-bottom depth: 45 + 114mm (1.77 + 4.49in)