

DEFIBRILLATOR ALARMED CABINET - 147SM-1

The Modern Metal alarmed AED wall cabinet is perfect for storing and maintaining the location of your AED.

The Aed Authority alarmed defibrillator cabinet is the perfect way to safely secure your defibrillator in any environment. Bright in colour and marked with AED signage this cabinet is easily distinguishable in an emergency situation. The AED Wall Cabinet suits almost all AED units on the market.

Dimensions are 400mmH x360mmW x 200mmD



KEY FEATURES:

- White powder coated cabinet
- High intensity red LED light on top of cabinet
- Loud 80dB alarm sounds and red light flashes when door is opened
- Marked with internationally recognised signage and 000 number
- Suits almost all Automated External Defibrillators
- Pre Drilled mounting holes at rear of cabinet

80
Decibel
Alarm tone