

ELECTRODE PRODUCT GUIDE



Product Name	Differentiation Features	Location Used	Defibrillator
CPR-D-padz 	<ul style="list-style-type: none"> Enables Real CPR Help and See-Thru CPR One-piece design simplifies proper positioning on the chest One size fits all, with ability to expand for larger chest 	<ul style="list-style-type: none"> EMS & Fire Public Access 	<ul style="list-style-type: none"> M Series® E Series® R Series® AED Plus® AED Pro® X Series™
CPR Stat-padz 	<ul style="list-style-type: none"> Enables Real CPR Help and See-Thru CPR Multi-function, two-piece design ensures all therapies can be delivered Accommodates AP and AA pad placement 	<ul style="list-style-type: none"> Hospital Med/surge floors ER & ICU EMS & Fire 	<ul style="list-style-type: none"> M Series E Series R Series AED Plus AED Pro X Series
OneStep Basic 	<ul style="list-style-type: none"> Multi-function electrode Automated self-test with R Series 	<ul style="list-style-type: none"> Hospital 	<ul style="list-style-type: none"> M Series R Series X Series
OneStep Pacing 	<ul style="list-style-type: none"> OneStep Basic features plus: OneStep Pacing 	<ul style="list-style-type: none"> Hospital 	<ul style="list-style-type: none"> M Series R Series X Series
OneStep CPR (AP) 	<ul style="list-style-type: none"> OneStep Basic features plus: Real CPR Help See-Thru CPR Anterior - Posterior placement OneStep Pacing 	<ul style="list-style-type: none"> Hospital 	<ul style="list-style-type: none"> M Series R Series X Series
One Step CPR (AA) 	<ul style="list-style-type: none"> OneStep Basic features plus: Real CPR Help See-Thru CPR Anterior - Anterior placement OneStep Pacing 	<ul style="list-style-type: none"> Hospital 	<ul style="list-style-type: none"> M Series R Series X Series
OneStep Complete 	<ul style="list-style-type: none"> All OneStep Basic, Pacing, and CPR capabilities 	<ul style="list-style-type: none"> Hospital 	<ul style="list-style-type: none"> M Series R Series X Series
OneStep Pediatric 	<ul style="list-style-type: none"> Designed for patients 0 - 8 years of age and less than 55 pounds (25 kg) Automatically adjusts joules setting to pediatric levels 	<ul style="list-style-type: none"> Hospital 	<ul style="list-style-type: none"> M Series R Series X Series
OneStep Pediatric CPR 	<ul style="list-style-type: none"> Designed for patients 0 - 8 years of age and less than 55 pounds (25 kg) Automatically adjusts joules setting to pediatric levels. Real CPR Help See-Thru CPR 	<ul style="list-style-type: none"> Hospital 	<ul style="list-style-type: none"> M Series R Series X Series
Pedi-padz® II 	<ul style="list-style-type: none"> Connects directly to AEDs & ALS devices Pediatric intelligence function with ZOLL AEDs 	<ul style="list-style-type: none"> EMS & Fire Public Access 	<ul style="list-style-type: none"> M Series E Series R Series AED Plus AED Pro X Series
Pedi-padz Multi-function Liquid or Solid Gel 	<ul style="list-style-type: none"> Multi-functional electrode for pediatric patients Available in liquid or solid gel 	<ul style="list-style-type: none"> Hospital EMS & Fire 	<ul style="list-style-type: none"> M Series E Series R Series X Series
Pedi-padz Radiolucent Solid Gel 	<ul style="list-style-type: none"> For use in Cath and RP labs for procedures Virtually invisible to X-ray 	<ul style="list-style-type: none"> Hospital Cardiac catheterization lab Electrophysiology labs 	<ul style="list-style-type: none"> M Series E Series R Series X Series
Pro-padz® Multi-function Liquid or Solid Gel 	<ul style="list-style-type: none"> Available in liquid or solid gel Available with or without tray Used in Cath and EP labs for procedures 	<ul style="list-style-type: none"> Hospital Cardiac catheterization lab Electrophysiology labs 	<ul style="list-style-type: none"> M Series E Series R Series X Series
Pro-padz Radiolucent Liquid or Solid Gel 	<ul style="list-style-type: none"> For use in Cath and EP labs for procedures - Virtually invisible to x-ray Pre-Connect (solid only) Speed pack 	<ul style="list-style-type: none"> Hospital Cardiac catheterization lab Electrophysiology labs 	<ul style="list-style-type: none"> M Series E Series R Series X Series
Pro-padz Sterile 	<ul style="list-style-type: none"> For use in operating rooms Helps maintain sterile barrier Packaged in two pouches Available in 54" or 10 feet lead wires 	<ul style="list-style-type: none"> Hospital Operating rooms 	<ul style="list-style-type: none"> M Series E Series R Series X Series
Stat-padz® 	<ul style="list-style-type: none"> Multi-function Pre-Connect Speed pack Different pad shapes and colors 	<ul style="list-style-type: none"> EMS & Fire Hospital 	<ul style="list-style-type: none"> M Series E Series R Series X Series
Stat-padz II 	<ul style="list-style-type: none"> Connects directly to AEDs and ALS devices 	<ul style="list-style-type: none"> EMS & Fire Hospital 	<ul style="list-style-type: none"> M Series E Series R Series AED Plus AED Pro X Series

See reverse side for proper cable/connector requirements.

CPR ELECTRODE GUIDE



CPR-D-padz®



CPR stat-padz®



OneStep™ Complete



OneStep™ CPR AP




OneStep™ CPR AA



OneStep™ Pediatric CPR



AED Plus®



- Connects directly for Real CPR Help® and defibrillation

- Connects directly for Real CPR Help and defibrillation

- Not compatible

- Not compatible

- Not compatible

- Not compatible

AED Pro®



- Connects directly for Real CPR Help, See-Thru CPR and defibrillation

- Connects directly for Real CPR Help, See-Thru CPR and defibrillation

- Not compatible

- Not compatible

- Not compatible

- Not compatible

X Series™



- Requires CPR Connector (PN# 8000-0370) for CPR Feedback
- Connects directly for defibrillation

- Requires CPR Connector (PN# 8000-0370) for CPR Feedback
- Connects directly for defibrillation

- Connects directly for all capabilities

- Connects directly for all capabilities

- Connects directly for all capabilities

- Connects directly for all capabilities

E Series®



- Requires CPR Connector (PN# 8000-0370) for CPR Feedback
- Connects directly for defibrillation

- Requires CPR Connector (PN# 8000-0370) for CPR Feedback
- Connects directly for defibrillation

- Not compatible

- Not compatible

- Not compatible

- Not compatible

R Series®



- Requires R Series® CPR Connector for CPR (PN# 8009-0020)
- Connects directly for defibrillation

- Requires R Series CPR Connector for CPR (PN# 8009-0020)
- Connects directly for defibrillation

- Connects directly for all capabilities

- Connects directly for all capabilities

- Connects directly for all capabilities

- Connects directly for all capabilities

Shelf Life

5 yrs

2 yrs

2 yrs

2 yrs

2 yrs

2 yrs

Quantity and Part Numbers

Individual: 8900-0800-01

Individual: 8900-0402
8 per case: 8900-0400

Individual: 8900-0224-01
8 per case: 8900-0214-01

Individual: 8900-0223-01
8 per case: 8900-0213-01

Individual: 8900-0225-01
8 per case: 8900-0217-01

Individual: 8900-0219-01
8 per case: 8900-0220-01