



HEARTSINE AED TRAINERS



REALISTIC TRAINING

The easy-to-use HeartSine samaritan PAD Trainer with remote control guides users through simulated analysis, energy delivery and prompted cardiopulmonary resuscitation (CPR) with the look and feel of a live HeartSine samaritan public access defibrillator (PAD) without the actual charge and discharge of an electrical shock.

The HeartSine samaritan PAD Trainer is available in three models: SAM 450P, SAM 360P and SAM 350P. The SAM 450P Trainer simulates the SAM 450P AED by providing the same visual and verbal CPR Rate Advisor prompts.

SIMULATE THE FULL RANGE OF “REAL LIFE” SCENARIOS

The multi-function remote control enables you to select among six AHA/ERC scenarios for live demonstrations or training. Additionally, the SAM 450P Trainer remote control provides the training instructor with the flexibility to choose the appropriate CPR Rate Advisor prompts based on the student's CPR compression rate.

Each Trainer is operated by a specific remote control, enabling multiple devices to be used in close proximity without interference. The remote control will operate with the Trainer within a range of 3 to 5 metres (10 to 15 feet).

ACCESSORIES



Remote for use with the HeartSine Samaritan Trainer. allows instructor to simulate various scenarios.

HeartSine samaritan PAD Trainer remote



HeartSine Samaritan Trainer, includes 1 replacement electrode cartridge.

Trainer-Pak



HeartSine Samaritan Trainer, includes 1 set of replacement Trainer Electrodes.

Replacement electrodes



Battery charger for use with the HeartSine Samaritan Trainer.


Battery charger



USB Data Cable for use with the HeartSine Samaritan Trainer.

USB data cable





03 8710 8666

sales@aedauthority.com.au

41 Marconi Drive, Dandenong South VIC 3175



aedauthority.com.au