

WIDE RANGE

ZOLL® offers a wide variety of accessories for defibrillation, monitoring, and multifunction use. Designed to accommodate various clinical needs, ZOLL's electrodes are easy to apply and use whilst the batteries that power ZOLL® devices are reliable and easy to maintain.



For use with the ZOLL AED Plus® and AED Pro® defibrillators, includes 1 set of adult electrodes and replacement instructions.

CPR-D-padz™



For use with the ZOLL AED 3 defibrillator, includes 1 set of adult electrodes and replacement instructions.

CPR Uni-padz™



For use with the ZOLL AED Pro® defibrillator, includes 1 set of adult electrodes and replacement instructions.

CPR Stat-padz® (paediatric electrode pads)



For use with the ZOLL AED Plus® and AED Pro® defibrillators, includes 1 set of adult electrodes and replacement instructions.

Pedi-padz® II





For use with the ZOLL AED Pro® defibrillator, includes 1 set of ECG electrodes and replacement instructions.

ECG Electrodes



10 pack of Duracell Lithium Batteries for the Zoll AED Plus® defibrillator.

Type 123 Lithium Batteries (Pack of 10)



For use only with the Zoll AED 3 defibrillator, includes 1 battery and replacement instructions.

ZOLL AED 3 battery



Non rechargeable battery for use only with the Zoll AED Pro® defibrillator, includes 1 battery and replacement instructions.

AED Pro Lithium-ion Pack (Non-rechargeable)



Rechargeable battery for use only with the Zoll AED Pro® defibrillator, includes 1 battery and replacement instructions.


SurePower® Battery Pack (Rechargeable)





For use only with the Zoll AED Pro® defibrillator rechargeable batteries, charges 1 battery.

SurePower® Single-Bay Battery Charger



For use only with the Zoll AED Pro® defibrillator rechargeable batteries, charges up to 4 batteries.

SurePower® Charging Station



ECG cable for the Zoll AED Pro® defibrillator.

AED Pro ECG Cable AAMI



Soft carry bag for the Zoll AED Plus® defibrillator.

AED Plus Black Carry Bag



Soft carry case for the Zoll AED 3 defibrillator.

ZOLL AED 3 carry case



SPECIFICATIONS

Zoll Batteries

Battery capacity:	Typical new battery running at an ambient temperature of +20° C to +25° C (68° F to 77° F) can provide: 140 defibrillator discharges at maximum energy (200 joules), or 6 hours of continuous monitoring (with 2-minute CPR periods) Note: CPR periods shorter than 2 minutes can decrease the operating time that can be obtained from a new battery
Type:	Disposable, sealed lithium manganese dioxide
Battery standby life (once installed):	5 years with weekly self-test. End of life indicated by blank status window
Battery shelf life:	Store for up to 2 years at 23°C (77°F) prior to installing to maintain above battery life
Temperature:	0°C to 50°C (32°F to 122°F)
Humidity:	10% to 95% (non-condensing)
Height:	27.75 mm (1.0 in)
Width:	133 mm (5.16 in)
Depth:	88 mm (3.5 in)
Weight:	317.5 grams (0.7 lbs)
Nominal voltage:	12 volts

CPR-D-padz™

Shelf life:	5 years
Conductive gel:	Polymer Hydrogel
Conductive element:	Tin
Packaging:	Multilayer foil laminate pouch
Impedance class:	Low
Cable length:	91cm (36 in)
Design standards:	Meets requirements of IEC 60601-2-4, 2010

CPR Stat-padz®


Shelf life:	2 years
Conductive gel:	Polymer Hydrogel
Conductive element:	Tin
Packaging:	Multilayer foil laminate pouch
Impedance class:	Low
Cable length:	124.5 cm (49 in)
Design standards:	Meets requirements of IEC 60601-2-4, 2010

CPR Uni-padz™

Shelf life:	5 years
Conductive gel:	Polymer Hydrogel
Conductive element:	Tin
Packaging:	Multilayer foil laminate pouch
Impedance class:	Low
Cable length:	142 (+/-3.8) cm; 56 (+/- 1.5) in
Design standards:	Meets requirements of IEC 60601-2-4, 2010

Pedi-padz® II

Age for use:	0-8 years
Weight for use:	Under 25 kg (55 lbs)
Shelf life:	5 years
Conductive gel:	Polymer Hydrogel
Conductive element:	Tin
Packaging:	Multilayer foil laminate pouch
Impedance class:	Low
Cable length:	86.4 cm / 34.0 inches
Design standards:	Meets requirements of IEC 60601-2-4, 2010



03 8710 8666

sales@aedauthority.com.au

41 Marconi Drive, Dandenong South VIC 3175



aedauthority.com.au