











INDOOR ALARMED DEFIBRILLATOR CABINET

The modern metal alarmed indoor AED wall cabinet is the perfect way to store, maintain and mark the location of your AED.

The AED Authority alarmed indoor defibrillator cabinet is the perfect way to safely secure your defibrillator in any indoor environment. Bright in colour and marked with internationally recognised signage this cabinet is easily distinguishable in an emergency situation.

The AED Wall Cabinet suits almost all AED units on the market.

KEY FEATURES:

-  White powder coated cabinet
-  **Height:** 400mm
-  **Width:** 360mm
-  **Depth:** 180mm
-  High intensity red LED light on top of cabinet
-  Loud 80DB alarm sounds and red light flashes when door is opened
-  Marked with internationally recognised signage and 000 number
-  Suits almost all Automated External Defibrillator models
-  Pre drilled holes for easy installation
-  Battery operated



03 8710 8666

sales@aedauthority.com.au

41 Marconi Drive, Dandenong South VIC 3175



aedauthority.com.au