

DEFIBRILLATOR COMPARISON



HeartSine Samaritan 360P



Lifepak CR2 Essential

RRP	\$2,495.00	\$3,095.00
Mode of Operation <i>Auto: Shock delivered automatically</i>	Auto	Auto
Carry Case	Soft case	Hard shell case
First Responders Kit <i>CPR Mask, Scissors, Razor, Hand wipes & Gloves</i>	Included	Included
ANZCOR <i>Meets Australian Resuscitation Council Guidelines</i>	No	Yes
Metronome <i>Audible guide to pace CPR</i>	Yes	Yes
CPR Coaching <i>Audible guidance</i>	Yes	Yes
Paediatric/Child Mode <i>Reduced energy delivery</i>	Pads available - Additional cost	Child button on unit
Shock Delivery/Joules <i>Rate of energy delivered during shock</i>	150 - 150 - 200	200 - 300 - 360
Warranty - Years	8	8
Battery & Electrode Life - Years	4	4
Battery & Electrode Replacment Cost	Adult Battery Pad-Pak - \$230.00 Paediatric Battery Pad-Pak - \$290.00 <i>(Integrated battery pad-pak)</i>	Battery - \$330.00 Electrode Pad - \$130.00
Weight	1.1kg	2.0kg
IP Rating <i>Dust & water protection</i>	IP56	IP55